

## The unique Vivatmo system allows for efficient asthma management

The Vivatmo product portfolio consists of the

- ▶ **Vivatmo** *pro* professional device for healthcare professionals
- ▶ **Vivatmo** *me* home-use device for asthma patients
- ▶ **Vivatmo** *app* digital asthma diary for documentation of patient's FeNO values, medication, and health condition



Therapy monitoring

Improved adherence

Better asthma management

#### **Benefits of the Vivatmo system**

- ▶ Optimum asthma monitoring in practices and at patient's home
- ▶ FeNO readings as valuable information for therapeutic decisions and therapy monitoring
- ▶ Appreciative cooperation between doctors and patients results in better adherence and overall satisfaction
- ► High-precision sensor technology and Bosch quality made in Germany

## What makes Vivatmo special

- **1. Perfect for routine use:** The Vivatmo devices are **maintenance- and calibration-free.** They are easy, quick, and intuitive to use for beginners and experts.
- **2. Takeaway:** The Vivatmo handheld is **mobile**, **wireless**, light-weight and particularly ergonomic.
- **3.** More than just a mouthpiece: The ergonomic design of the disposable mouthpiece supports a correct breathing maneuver. The breath sample is filtered, cleaned, and prepared for NO measurement. Since users only exhale into the mouthpiece, cross-infection risk is greatly reduced.

#### Vivatmo pro for healthcare professionals

# FeNO measurement for higher throughput in practices and clinics

- ▶ Up to 5,000 measurement trials or 5 years lifetime without any maintenance or calibration
- ► Base station with large high-quality touchscreen and mobile handheld for stand-alone measurements
- ▶ Different kinds of **visualization** designed to support patient's measurement performance as optimally as possible
- Quality control with two reference points and ambient NO measurement



# Easy and secure data management and exchange

- ► Patient data management system of up to 10,000 patient files, including export function and immediate printing
- ► Backup functionality and optional user administration with Login to support data security in clinical environment
- ► Technical support of several modern data exchange options: Wireless using Wifi, via USB or Ethernet ports, using the clinic's EPR via HL7/ GDT
- ► Remote software updates available via Vivasuite connectivity platform



## Vivatmo me for asthma patients

## The world's first FeNO measuring device for home use

- ► Mobile handheld device combined with Vivatmo app for better asthma management
- ► Up to 1,000 measurements or 3 years lifetime without any maintenance or calibration
- ► Easy and intuitive user guidance during the measurement with three-color LED and feedback sound
- ► **Result interpretation** is supported by traffic light classification and trend arrows

# Ideal supplement: free Vivatmo app as digital asthma diary

► The app gives a good overview of FeNO values, peak flow results, medications, and asthma attacks

► Regional pollen forecast is provided in cooperation with the University of Vienna

► Monthly reports of all data can be saved or shared with doctor showing the course of asthma disease



#### **Bosch Healthcare Solutions GmbH**

Stuttgarter Str. 130 D-71332 Waiblingen Tel.: +49 -(0) 711-811-0

www.vivatmo.com

#### Distributed by:



#### COSMED Srl

Via dei Piani di Monte Savello 37 Albano Laziale - Rome 00041 Italy +39 (06) 931-5492 Phone +39 (06) 931-4580 Fax

cosmed.com

Technical data subject to modification without notice.

Printed in Germany. Any color variations are due to the limitations of the printing process.

Printed on 100% chlorine-free bleached cellulose manufactured in an environmentally friendly process.

Not all products are available in all regions. Check with your local representative for availability in specific markets. For In Vitro Diagnostic Use Only.