



MONARK 871E



MAIN FEATURES

- Arm and leg ergometer
- Compact size

THE BASICS

- Good cycle feeling
- User friendly display
- Can be used from wheelchair
- No power adaptor needed
- Pedals and handles included

Medical Device class IIa

Specifications

MONARK 871E

Certifications	EU directive 93/42/EU class IIa
Frame	
Type	Arm ergometer
Handlebar	n/a
Stem	n/a
Seatpost	n/a
Saddle	n/a
BB	68/115 mm
Crank	Steel, 152 mm
Pedals	9/16", handles
RESISTANCE	
System	Monark Basic
Flywheel	12 kg
DISPLAY	
RPM	Yes
HR	n/a
Time	Yes
Speed	n/a
Distance	n/a
Watt	n/a
Connectivity	
RS232	n/a
Usb	n/a
WEIGHT & MEASUREMENTS	
Max user weight	n/a
Length [mm]	600 mm
Width [mm]	180 mm at front support leg, 750 mm at rear support leg
Height [mm]	530 mm
Weight [kg]	Ergometer 25 kg
POWER REQUIREMENTS	Batteries in the display



WHEELCHAIR

Adapted for use from a wheelchair

MEDICAL CE

Monark Exercise AB is since 1992 ISO 9001 certified. All models are:

- Approved under CE Certificate No. 15 36 01 (MDD IIa)
- Anti-corrosion treated and powder coated
- Made in Sweden

