Ŀ



MONARK 881E

MAIN FEATURES

- Arm and leg ergometer
- Individually adjustable crankarms
- Possible to calibrate

THE BASICS

SPORTS & MEDICA

Good cycle feeling

MONARK

11

- User friendly display
- Can be used from wheelchair
- No power adapter needed
- Pedals and handles included

MONARK REHAB TABLE

- Adjustable hight
- Adapted for 881E
- Lockable
- Wheelchair size

Medical Device class IIa



Specifications MONARK 881E

Certifications	FU directive
	93/42/EU class IIa
Frame	
Туре	Arm ergometer
Handlebar	n/a
Stem	n/a
Seatpost	n/a
Saddle	n/a
BB	n/a
Crank	Steel, 152 mm
Pedals	1/2"
RESISTANCE	
System	Spring
Flywheel	11 kg
DISPLAY	
RPM	Yes
HR	n/a
Time	Yes
Speed	n/a
Distance	n/a
Watt	n/a
Connectivity	
RS232	n/a
Usb	n/a
WEIGHT & MEASUREMENTS	
Max user weight	n/a
Lenght (mm)	540 mm
Width (mm)	470 mm
Height (mm)	550 mm
Weight (kg)	22 kg
POWER REQUIREMENTS	Batteries in the display



WHEELCHAIR Adapted for use from a wheelchair

MEDICAL CE

Monark Exercise AB is since 1992 ISO 9001 certified. All models are:

- Approved under CE Certificate No. 153 601 (MDD IIa)
- Anti-corrosion treated and powder coated
- Made in Sweden

Upplaga 1. 2016-02-28



