



MONARK 927E



MAIN FEATURES

- User friendly exercise bike for rehabilitation
- Easy entry frame
- Extra durable design

THE BASICS

- Good cycle feeling
- 5 kHz HR systems
- User friendly display
- No power adaptor needed

Max user weight 250 kg

Medical Device class IIa

Specifications

MONARK 927E

Certifications	EU directive 93/42/EU class IIa
Frame	
Type	Medical Frame
Handlebar	Monark Handlebar
Stem	n/a
Seatpost	Monark Standard
Saddle	Saddle, black
BB	68/122 mm
Crank	Steel, 170 mm
Pedals	9/16", pedals with foot straps
RESISTANCE	
System	Monark Basic
Flywheel	12 kg
DISPLAY	
RPM	n/a
HR	5 kHz
Time	Yes
Speed	Yes
Distance	Yes
Watt	n/a
Connectivity	
RS232	n/a
Usb	n/a
WEIGHT & MEASUREMENTS	
Max user weight	250 kg
Length [mm]	1240 mm
Width [mm]	500 mm at handlebar, 640 mm at support legs
Height [mm]	1260 mm at handlebar 780-1170 mm at saddle
Weight [kg]	46 kg
POWER REQUIREMENTS	Batteries in the display



HEART RATE

Equipped to measure heartrate via chestbelt.

MEDICAL CE

Monark Exercise AB is since 1992 ISO 9001 certified. All models are:

- Approved under CE Certificate No. 15 36 01 (MDD IIa)
- Anti-corrosion treated and powder coated
- Made in Sweden

