



MONARK 939 NOVO

MAIN FEATURES

- Designed for external control
- Monark novo touch display
- Monark novo Control unit
- Super fast motor
- Easy entry frame

THE BASICS

- Speed independent and dependent power control
- Easy to control and calibrate
- Can be controlled by two external devices (RS232 or USB) or bike display
- Resistive LCD Touch screen
- BLE Smart and ANT+ HR systems
- Easy to connect to external devices
- Workload range 4-1400 W (200 rpm)
- Protocols: Depending external device functionality
- Training mode: Watt, KPM/min, VO2, kP and N



Max user weight 250 kg

Medical Device class IIa

Monark Test Software

