## **SPORTS** RANGE

# MONARK LC7 TT novo

#### MAIN FEATURES

- Designed for TimeTrial/Open-end testing
- Workload control from bars or PC/ECG
- Monark novo touch display
- Monark novo control unit
- Super fast motor
- Perfect fit frame for true cycling position

#### THE BASICS

- Perfect cycle feeling
- Easy to use your own saddle, handlebar and pedals
- Speed independent and dependent power control
- Handlebars are easy to change to suit different tests
- Easy to control and calibrate
- Can be controlled by two external devices (RS232 or USB)
- BLE Smart and ANT+ HR systems
- Easy to connect to external devices
- Workload range 4-1400 W (200 rpm) (8-2800 W option)
- Protocols: Full Monark Test Software functionality
- Traning mode novo display: Watt, KPM/min, VO2, kP and N

Max user weight 180 kg Medical Device class IIa Mo

Monark Test Software



nNn

# Specifications MONARK LC7 TT novo

Certifications	EU directive 93/42/EU class IIa
Frame	
Туре	Sports, perfect fit
Handlebar	Racing handlebar
Stem	Adjustable, horizontal 200 mm, vertical 500-910 mm
Seatpost	Adjustable, horizontal 200 mm, vertical 530-940 mm
Saddle	Racing
BB	68/107 mm
Crank	Steel, 172.5 mm
Pedals	9/16", combi SPD
RESISTANCE	
System	Pendulum
Flywheel	20 kg
DISPLAY	
RPM	Yes
HR	ANT+, bluetooth smart
Time	Yes
Speed	Yes
Distance	Yes
Watt	Yes
Connectivity	
RS232	Yes
Usb	Yes
WEIGHT <del>8</del> MEASUREMENTS	
Max user weight	180 kg
Length (mm)	1405-1630 mm
Width (mm)	640 mm
Height (mm)	1500 mm at display /PC 1240 mm at saddle
Weight (kg)	77 kg
POWER REQUIREMENTS	220V (12V, 5,0A)



#### ECG

Can be controlled from different ECG brands either analogue or digital.



#### PC CONNECTION

Direct connection to PC through a serial cable.



#### SOFTWARE

For fitness testing and in some cases, even training can be used with this model.



**HEART RATE** Equipped to measure heartrate via chestbelt.



#### TEST

Offers the possibility of performing different tests thanks to the accurate workload which can be calibrated.



### CONSTANT POWER MODE





FIXED PROTOCOL

For test or training controlled through the display or PC program.

#### MEDICAL CE

Monark Exercise AB is since 1992 ISO 9001 certified. All models are:

- Approved under CE Certificate No. 15 36 01 (MDD IIa)
- Anti-corrosion treated and powder coated
- Made in Sweden





